

FT-IR Investigation of Glucoheptonate Chelation


AAPFCO Winter 2020 Meeting, New Orleans LA

Dr. Eric Johnson
Chemical Dynamics, Inc.
4206 Business Lane
Plant City FL 33566

CHEMICAL DYNAMICS, INC.

LIGHT


Classification of Electromagnetic Radiation



Light rays are electromagnetic waves: they can be define by their frequency or their wavelength (λ).


Infrared Spectroscopy: Overview

- Chemical bonds are like springs. They vibrate (bend, stretch, etc).
- Infrared light can be absorbed by these vibrations, increasing their energy levels.
- Frequency is characteristic of the bonds involved
- polarity, geometry, symmetry, selection rules
- Spectral range typically 4000 to 450 cm^{-1} .
- (2500 nm to 22200 nm)
- Molecular fingerprint



Attenuated Total Reflection - ATR

- When radiation is totally reflected internally at an interface with a material of lower refractive index the electric field penetrates a short distance beyond the interface. This field, called an *evanescent wave*, can be absorbed by a sample on the surface.
- The electric field decreases away from the surface and the *penetration depth* (d_p) is defined as the distance at which it falls to $1/e$ of its value at the surface.
- Penetration depth depends on the angle of incidence (θ) and the refractive indices of crystal and sample and is also proportional to the wavelength (1).
- For a diamond ATR element, the penetration depth at 1000 cm^{-1} is typically $1\text{--}2\text{ }\mu\text{m}$.
- Because of the wavelength dependence, the penetration depth increases by a factor of 10 between 4000 and 400 cm^{-1}



[1] F.M. Mirabella, *Internal Reflection Spectroscopy* (Marcel Dekker, New York, 1993).

Perkin Elmer Spectrum 2 with UATR



- Single reflection diamond ATR crystal, epoxy bonded into 316 Stainless Steel
- Minimal Sample prep, just evaporate a drop of solution onto crystal
- Easy clean up
- A few minutes per sample

Studies:

- Prepare solutions of metal cations and sodium glucoheptonate.
- Zn, Mn, Fe and Cu
- Sodium Glucoheptonate is not chelated² while Zinc Glucoheptonate, Zinc Gluconate^{3,4} and Manganese Gluconate⁵ are.
- Compare FT-IR ATR spectra of glucoheptonates of Mn, Cu and Fe with zinc glucoheptonates.

2. 2. Y.J Park, B.H Lee, *Acta Cryst.* C57, 12 (2001).

3. H.A. Tajmir-Riahi, *Can. J. Chem*, **67**, 651 (1988)




4. G.M. Escandar, M.G. Sierra, J.M.S. Peregrin, G. Labadie, M. Santoro, A. Frutos, L.F. Sala, *Polyhedron*, **13**, 909 (1994).

5. T. Lis, *Acta Cryst.*, **B35**, 1699 (1969).

Experimental

- Prepare series of aqueous solutions with constant sodium glucoheptonate concentration while varying cation content.
- Cations derived from metal chlorides. Byproduct, sodium chloride, does not have covalent bonds and is infrared transparent.
- Solution is placed on ATR crystal (after background of clean crystal obtained) and dry air used to evaporate water leaving behind metal glucoheptonate.
- Spectra recorded in absorbance and normalized with respect to glucoheptonate content.
- Grams/AI (Thermo) used for analysis

What Kind of Changes to be seen?



- Reacting glucoheptanoic acid to the sodium salt: Conversion of carboxylic acid to carboxylate ion, chelation not present

Chelation with metal cation:


- Change in frequency and intensity of carboxylate ion⁶
- Change in C-O bending / deformation modes due to change in bond angle (5 membered ring)^{7,8}
- Changes in O-H stretching frequency and Intensity^{7,8}

6. H.A. Tajmir-Raihi, *Carbohydr. Res.*, **122**, 241(1983)

7. E. Wiercigroch, E. Szafraniec, K. Czamara, M. Z. Pacia, K. Majzner, K. Kochan, A. Kaczor, M. Baranska, Kamilla Malek, *Spectrochimica Acta Part A: Molecular and Biomolecular Spectroscopy*, **185**, 317 (2017)


8. A.P. Kilimov, M.A. Svechnikova, V.I. Shevchenko, V.V. Smirnov, F.V. Kvasnyuk-Mudri, S.B. Zotov, *Chemistry of Heterocyclic Compounds*, **3**, 647 (1967).

What Kind of Changes to be seen?




Absorbance / Wavenumber (cm-1)

Control: Zinc Glucoheptonate




Manganese Glucoheptonate




Copper Glucoheptonate

Complex mixture of multiple types of chelated species.



9. G.M. Escandar, L.F. Sala., "Complexes of Cu(II) with D-almonic and D-alduronic acids in aqueous solution" *Can. J Chem.*, **70**, 2053 (1992)

Copper Glucoheptonate




Absorbance / Wavenumber (cm-1)

2/20/2020

Chemical Dynamics, Inc.

13

Iron Glucoheptonate



Summary

- ATR effective tool for investigating changes related to metal glucoheptonate chelation.
- Zinc Glucoheptonate model compound
- Manganese Glucoheptonate follows similar spectral changes
- Copper complex mixture of chelated forms
- Iron similar to copper

Limitations and Opportunities for More Work

- Work is primarily qualitative.
- Compare sodium glucoheptonate with metals to reaction with glucoheptanoic acid?
- Chemical environment effects on spectra?

Acknowledgements

- Ray Thompson, Perkin Elmer
- Kolin Alfonso

Questions?